

PHERO-stab 0305 E

The new power supply in the approved PHERO-stab series !



Technical data:

Voltage (fixed):	50 / 70 / 80 / 90 / 100 / 110 / 120 125 / 175 / 200 / 225 / 250 / 275 / 300
Current:	max. 50 mA
Timer:	1....99 Min. or permanent
Ports:	2 Parallel
Dimensions:	20 x 9 x 15 cm
Order-code:	131112

Specials:

- The PHERO-stab 0305 E works with constant voltage. If the maximum current from 50 mA is reached, the power supply changes to constant current. This is shown with LED's.
- It is possible to work with or without the timer.
- After the adjusted time, 20% of the selected voltage rested, till it is switched off. In this way the samples will not diffuse.